


## The role of sorting nexin-1 (SNX1) during *Salmonella* infection




**Figure 1:** The boxed region in part A shows that SNX1 (green) accumulates at membrane ruffles (Scale bar: 10  $\mu\text{m}$ ). A magnification of this boxed region is shown in part B, in which different orientations of the 3D rendered view are displayed, with the green (SNX1) channel omitted in B(iii). Part C shows a single optical z-section of the maximum projections shown in figures A and B, with insets providing the respective YZ- and XZ-views.

Professor Peter Cullen and colleagues are interested in the role of the mammalian sorting nexins in host-pathogen interactions. *Salmonella enterica* serovar Typhimurium (*S. Typhimurium*) is the leading cause of human gastroenteritis and invades non-phagocytic cells by actively inducing its uptake. Once inside, *S. Typhimurium* establishes a replication niche within membranous *Salmonella*-containing vacuoles (SCVs) which are enriched in phosphatidylinositol 3-monophosphate [PtdIns(3)*P*] and 'blend in' with the host's endosomal system.

In this work the researchers use a combination of live cell and immunofluorescence microscopy (using the PerkinElmer UltraVIEW LCI confocal Optical Scanner) to show that sorting nexin-1 (SNX1) is crucially involved in remodeling nascent SCVs. SNX1 contains a sorting nexin phox (PX)-homology domain that binds to PtdIns(3)*P* and a membrane-binding module (called a BAR domain) which induces membrane tubulation.

MDCK cells were infected with *S. Typhimurium* strain SL1344 for 15 minutes, then fixed, immunolabeled with anti-SNX1 (green), and treated with TRITC-phalloidin (red) and DAPI (blue). Ten optical z-slices were acquired ( $z=5.37 \mu\text{m}$ ) and data sets were deconvolved, 3D rendered and animated using Volocity software (Figure 1; see also <http://jcs.biologists.org/content/vol121/issue12/images/data/2027/DC1/018432-Movie1.mov>). SNX1 was shown to be readily recruited to sites of *Salmonella* invasion, and was particularly found clustered underneath actin-rich membrane ruffles. This recruitment of SNX1 was found to induce the formation of numerous long-range tubules. In order to better visualize the full extent of tubulation, cells expressing GFP-SNX1 were infected with SL1344, fixed and treated with TRITC-phalloidin and DAPI. Confocal imaging was used to acquire 29 optical z-slices ( $z=4.76 \mu\text{m}$ ) and data sets were deconvolved, 3D rendered, animated and exported using Volocity software (Figure 2; see also <http://jcs.biologists.org/content/vol121/issue12/images/data/2027/DC1/018432-Movie3.mov>). The results revealed an extensive tubular meshwork of GFP-SNX1 around the invading bacteria 15 minutes after infection.

The authors went on to show that SNX1-induced tubulation is required for SCV biogenesis as suppression of SNX1 leads to defects in the intracellular progression and replication of the invaded *Salmonella*.



**Figure 2:** MDCK cells containing a GFP-SNX1 fusion protein (green) were fixed after a 15 minute infection with SL1344, and treated with TRITC-phalloidin (red) and DAPI (blue). Images of a deconvolved stack are shown. GFP-SNX1 creates a meshwork around bacteria.