

Topic

This technical note outlines the Windows and Macintosh system requirements for the Improvisation License Server.

Discussion

Recommended Specification:

<i>Processor</i>	Macintosh G3, G4, G5 or Intel	Any Intel Pentium II CPU or higher. Any AMD Athlon CPU or higher Note: Intel Pentium or AMD K6/K6-2 CPUs of any speed are not supported.
<i>Operating System</i>	Mac OS 10.4.7 or higher	Windows Vista, Windows XP, Windows Server 2003 or Windows 2000 (Windows 95/98/ME or Windows NT 3.1/3.5/4.0 are not supported)
<i>RAM</i>	Apple recommends a minimum of 128 MB RAM to run Mac OS X	32 MB minimum, 128 MB preferred
<i>USB Ports</i>	One free USB port required for installation of dongle. A USB expansion card may be purchased for computers without built in USB.	
<i>Free hard disc space</i>	50 MB minimum, 100 MB preferred	

The server computer must be set to 'never sleep' and must be unlikely to be switched off.

The server computer may also act as a client running one of the applications. If the server computer is to also act as a client, the specification must be as for the client application. For clients to run Volocity see TN263: Volocity System Requirements.