

### Topic

This technical note outlines the Windows and Macintosh system requirements for the Improvision License Server.

### Discussion

#### Recommended Specification:

|                             |   |  |
|-----------------------------|---|--|
| <i>Processor</i>            | Macintosh G3, G4, G5 or Intel   | Any Intel Pentium II CPU or higher.<br>Any AMD Athlon CPU or higher<br>Note: Intel Pentium or AMD K6/K6-2 CPUs of any speed are not supported. |
| <i>Operating System</i>     | <b>Mac OS 10.4.7 or higher</b>  | <b>Windows Vista, Windows XP, Windows Server 2003 or Windows 2000</b><br>(Windows 95/98/ME or Windows NT 3.1/3.5/4.0 are not supported)        |
| <i>RAM</i>                  | Apple recommends a minimum of 128 MB RAM to run Mac OS X  | 32 MB minimum, 128 MB preferred  |
| <i>USB Ports</i>            | One free USB port required for installation of dongle.<br>A USB expansion card may be purchased for computers without built in USB. |  |
| <i>Free hard disc space</i> | 50 MB minimum, 100 MB preferred   |  |

The server computer must be set to 'never sleep' and must be unlikely to be switched off.

The server computer may also act as a client running one of the applications. If the server computer is to also act as a client, the specification must be as for the client application. For clients to run Volocity see TN263: Volocity System Requirements.